Mouth-actuated Tester GT-M03

Product Presentation



Application

This tester is suitable for the test of the safety of the Mouth-Actuated Toys, such as whistles with a ball, the sounders with a reed and other toys. This tester will blow-off certain volume gas in a certain speed and inhale the toy, blow-off and inhale alternately for 10 times, after the test, you should test the toy and the deciduous parts whether they are findings.

Standards

EN-71-1 section 8.17.2; ISO 8124-1 section 5.20 ASTM F963 section 8.13; GB 6675.2 section 5.20

Key Specification

Model	GT-M03
Blow and suck air volume	295±10 cm2
Blow and suck frequency	1 cycle/5s
Test air pressure	13.8 Kpa
Testing times	10 times
Air pressure meter	Digital, Units, Kpa; The positive and negative pressure display
Dimensions (L x W x H)	50 x 42 x 26cm
Weight	26Kg
Power supplier	AC220V 50/60HZ
Accessories	Tube gripper -3pcs



GESTER INTERNATIONAL CO.,LTD 7/F, No.15, Chifeng Road, Licheng Region, Quanzhou city of Fujian PR, 362000 China.

TEL: 86 0595-22515230 FAX: 86 0595 22515221

E-mail: sales3@gester-instruments.com Web: www.labtestingmachines.com